

# Drionic Iontophoresis for Hyperhidrosis

- Stop excessive sweating
- Units available for hands, feet and axillae
- Device uses tap water
- Use for 1/2 hour daily
- Usually stops sweating in 2 weeks
- Reusable
- Battery operated - use anywhere
- Replaceable batteries and pads
- Ergonomic construction



Tel: 071 138 8605  
Fax: 086 766 4462

Email: [sales@skintech.co.za](mailto:sales@skintech.co.za)  
[www.skintech.co.za](http://www.skintech.co.za)